

5TH EDITION.

PATENT
ABDOMINAL
and
OBSTETRIC

DOMEN

BELTS.
BELT CORSETS.
STOOP CURE
etc. etc.

E. R. WILLS,

80, PARK STREET,

BRISTOL.

WELLCOME INSTITUTE LIBRARY	
Coll.	welMOMec
Coll.	pam
No.	WE 158
	191 *
	W 74 P



22500837262

CONTENTS.

	PAGE
TESTIMONIALS	2
INTRODUCTION	3
SPECIAL NOTICE	4
SUPPORT WITHOUT PRESSURE.	5
“DOMEN” ABDOMINAL BELT SUPPORTS:—	
PREGNANCY BELT System I.	6
ACCOUCHEMENT BELT „ II.	7
HYPOGASTRIC BELT „ III.	8
BELT WITH SPINAL ATTACHMENTS „ IV.	9
UMBILICAL HERNIA BELT „ V.	10
FLOATING KIDNEY BELT „ VI.	11
INSTRUCTIONS FOR MEASURING	12
INSTRUCTIONS FOR WEARING	13
“DOMEN” ABDOMINAL BELT BANDAGES	14-16
“DOMEN” BELT CORSET	17
“DOMEN” STOOP CURE	18, 19
“DOMEN” MEN’S BODY BELTS	20
“DOMEN” CHOLERA BELTS	21
“DOMEN” ELASTIC HOSIERY, STOCKINGS, &c.	22, 23
“DOMEN” OCCASIONAL BELT	24

TESTIMONIALS, ETC.

A LARGE number of unsolicited Testimonials, to the "Domen" manufactures have been received from private persons and from members of the Medical Profession, referring in the highest terms to the efficiency and durability of the "Domen" Goods, and to the comfort derived from their wear.

A number of the Leading Hospitals, among which are
The London Hospital,
The Westminster Hospital,
The University College Hospital,
The South Devon and East Cornwall Hospital,
also use the "DOMEN" Specialities, with most beneficial results and to their entire satisfaction.

A FEW WORDS OF INTRODUCTION.

THE subject of Belts has lately been prominently before the Public, in the shape of so-called "Electric Belts," which, it is alleged, do not electrify. Even if some electric current had been generated, it is doubtful how far electrical influence can provide a remedy for the complaints which so-called "Electric Belts" have been advertised to cure. A succession of electric shocks may stimulate a nerve or a muscle which has through some cause become paralysed or disabled; but probably the useful purposes which electricity can serve are restricted to very narrow limits, and any cures effected by so-called "Electric Belts" have been mainly "faith cures."

The "DOMEN" BELTS COMPANY venture to recommend their perfected inventions to the Medical Profession and the Public on the grounds of careful construction, thorough adaptability to the most varying needs and requirements, and generally of honest excellence, free from all quackery and charlatanism.

Wholesale and Shipping Warehouse,
61, MOOR LANE, LONDON, E.C.

SPECIAL NOTICE.

EVERY "DOMEN" ARTICLE IS MADE OF THE BEST OBTAINABLE MATERIAL AND IN THE MOST APPROVED STYLE.

THE COMPANY HAVE MADE ARRANGEMENTS WITH DR. JAEGER'S SANITARY WOOLLEN SYSTEM COMPANY, LIMITED, LONDON, FOR KEEPING ALL "DOMEN" GOODS ON THEIR WELL-KNOWN AND HIGHLY RECOMMENDED SANITARY WOOLLEN SYSTEM.

WHEREVER THE JAEGER CLOTHING HAS BEEN ADOPTED, A SANITARY WOOLLEN BELT OR APPLIANCE IS LIKELY TO BE PREFERRED, BOTH FROM THE POINT OF VIEW OF HEALTH AND OF COMFORT.

THE "DOMEN" BELTS ARE GUARANTEED TO PERFECTLY FIT EVERY FIGURE IF THE MEASUREMENT BE CORRECTLY GIVEN ACCORDING TO PAGE 12.

THE "DOMEN" BELTS WILL NOT **Ruck Up** or **become Displaced** WHEN ONCE THEY HAVE BEEN PROPERLY PUT ON ACCORDING TO THE SIMPLE AND PLAIN INSTRUCTIONS PRINTED ON PAGE 13.

SUPPORT WITHOUT PRESSURE,

FOR

PREGNANCY, ACCOUCHEMENT, CORPULENCY.

THE "DOMEN" ABDOMINAL BELTS are a scientific improvement on the old-standing idea that the same belt will do for a number of greatly varying conditions. The whole subject of Abdominal Belts has now been carefully thought out, and every improvement that thorough anatomical knowledge and mechanical science can suggest has been adopted. Each special want or complaint has been provided with the appliance peculiarly suited to it, with the happiest and most beneficial results.


The "DOMEN" ABDOMINAL BELTS fit well, sit comfortably, and *do not get out of place*. They are made in six different forms, giving great comfort and relief in abdominal complaints, lessening the burden of pregnancy and the risk of miscarriage, and promoting a safe and easy delivery.

Such skilled and minute consideration of the remedies for different complaints has never before been placed at the disposal of Medical Men and their Patients in this country; and it is only necessary to examine the contents of this Catalogue to recognise the vast number of cases in which a Security, Relief, and Comfort hitherto unattainable can now be **immediately** and **cheaply** afforded.

6 Support without Pressure, "DOMEN" ABDOMINAL BELTS.

PREGNANCY BELT,

NOT TO BE WORN AFTER CONFINEMENT.



SYSTEM I.

Renders very effective service in Pregnancy, by greatly diminishing the suffering and burden usually entailed; much reduces the dangers attending delivery; the best preventive of premature confinement; affords instant relief in derangements of the womb, especially in prolapsus (falling of the womb), anteversion, and anteflexion. *The abdomen is literally carried*, and the muscles relieved from all excessive strain by the non-yielding material at the hips, and by the strong, broad band, which is kept in position by hooks as shown, and supports the lowest point of the body. The elastic front allows easy expansion where it is needed *and prevents undue pressure*. Prices for the various qualities, in all sizes up to 40 inches, are:—

A	W	E	K	J	R
Lavender Grey Jean.	Superior White Coutile	White Horsehair (Summer).	Black "Diamond Drill" Specialty.	On Dr. Jaeger's System. Natural Colour Lasting.	Special Horsehair Material.
List Nos. 701	702	703	704	705	706

"Obstetric Binders," for use immediately after delivery,
List Nos. 774 and 775.


NOTE.—Instructions for measuring and extra charges for sizes over 40 ins.
are given on page 12.

ACCOUCHEMENT BELT.

RESTORES FORMER SHAPELINESS

Pre-eminently useful after confinement, affording the most effectual support to the internal organs, preventing pressure on the fundus of the womb, and protecting against the numerous diseases and deformities which may occur from want of due precaution at these critical stages. The elastic sides of this Belt, combined with the non-yielding front, reduce the abdomen to its original size and shape-line; convalescence is thereby greatly accelerated, and weakness of the abdominal walls prevented.

SYSTEM II.



NOTE.—Although this Belt is not so eminently suitable as System I. for wear during pregnancy, it will do good service for those who are obliged to study economy.

Prices for the various qualities, in all sizes up to 40 inches, are:—

A	W	E	K	J	R	H
Lavender Grey Jean.	Superior White Coutile.	White Horsehair.	Black "Diamond Drill" Specialty.	On Dr. Jaeger's System, Natural Colour Lasting.	Special Horsehair Material.	Scarlet Wool Satin.
List Nos. 707	708	709	710	711	712	713

In A Quality this System is stocked in all sizes from 26 inches to 50 inches in circumference.

FIG. A.

AFTER OPERATION OF LAPAROTOMY.



After operation of Laparotomy this is the most efficient Belt that can be used, giving, when fitted with the velvet-covered pad, Fig. A, firm support to the weakened wall over and around the cicatrix. Isosceles triangularly shaped pads, slightly concave, velvet covered, of following sizes:— $5\frac{1}{2}$ by $3\frac{1}{4}$ ins., and $6\frac{1}{2}$ by 3 ins., and 7 by $4\frac{1}{2}$ ins. List No. 776. The pad is sewn fast in the required position on the inner side of the Belt.

System II. is also particularly serviceable in cases of corpulence or weakness of the abdominal regions for both **LADIES** and **GENTLEMEN**; its beneficial action will be most felt in *walking, riding, or during any active exercise.*

NOTE.—Instructions for measuring and extra charges for sizes over 40 inches are given on page 12.

HYPOGASTRIC BELT.

FIG. 1.




FIG. 2.



SYSTEM III.

This Belt should be employed in cases of extreme corpulency and of pendulous abdomen, when it will be found unsurpassed in point of efficiency and comfort. Its action is most beneficial when the womb is pressed up against the bladder, as also in the pains which so frequently make themselves felt in the hips, loins, and sides.

A peculiarly constructed flexible steel spring (Fig. 2), with a well-cushioned pad, maintains an equable and gentle pressure immediately above the pubes.

Prices, in all sizes up to 40 inches, are:—

Qualities kept in Stock.	A	E	J
	Lavender Grey Jean.	White Horsehair (Summer).	On Dr. Jaeger's System. Natural Lasting.
	List Nos. 714	715	716

In A Quality this System is stocked from 28 inches to 46 inches.

NOTE.—Instructions for measuring and extra charges for sizes over 40 inches are given on page 12.

ABDOMINAL SUPPORT WITH SPINAL SPRINGS.

FIG. 1.





FIG. 2.



SYSTEM IV. List No. 717.

The special construction of this Support is peculiarly adapted for relieving pains in the back and disorders caused by weakness of the spine.

Attached to the ends of this Support are two covered flexible steel springs (Fig. 2), the ends of which are furnished with well-cushioned pads, which act as comfortable and powerful stays to the spine.


In the case of an obese person with pendulous abdomen, the padded springs prevent the Support from cutting into the body and from pressing upon the spine.

Spinal springs, as illustrated in Fig. 2, can be attached to any of the Belts described on pages 5 to 10, and can be had separately. List No. 777.

System IV. Belts are supplied in all qualities in which System II. Belts are kept (see page 7).

NOTE.—Instructions for measuring and extra charges for sizes over 40 inches are given on page 12.

UMBILICAL HERNIA.



SYSTEM V.

This Support (Fig. 1) is the best appliance for Umbilical Hernia. In the case of an obese person with pendulous abdomen suffering from Umbilical Hernia, and especially with women past middle age, who have borne large families, the abdomen, from its relaxed condition, has a tendency to push an Umbilical Spring Truss out of position. In such instances the excellent shape and construction of this Belt maintain it in position, affording an exceedingly comfortable stay to the abdomen and the internal abdominal organs, without requiring it to be drawn too tight.

All Belts of this System are fitted with velvet-covered oval-shaped concave pads (Fig 2).

A similar pad with apex centre (Fig. 3) may be had instead. This should, however, only be used under medical authority.

Prices, in all sizes up to 40 inches, are :—

Qualities kept in Stock.	A Lavender Grey Jean.	E White Horsehair (Summer).	J On Dr. Jaeger's System Natural Lasting.
	List Nos 718	719	720

Stocked in sizes 28 inches to 50 inches in A Quality.

A special Belt for the above complaint has been designed for Men's use, which, however, is not intended for support to the abdomen, but merely for the hernia. Shape as Belt design, page 21, to buckle at back, with velvet-covered oval-shaped concave pad in centre. Made to order only. List No. 773.

NOTE.—Instructions for measuring and extra charges for sizes over 40 inches are given on page 12.

Support without Pressure, "DOMEN" ABDOMINAL BELTS. 11

FOR FLOATING KIDNEY.




FIG. 1.



FIG. 2.

SYSTEM VI.

Especially constructed for Displacement of the Kidneys, which it retains in their natural position, without injury or troublesome pressure, and with immediate cessation of pain. In ordering this Support, information should be given as to whether the displacement is on the right or the left, or on both sides. For the latter an extra charge of 2s. will be made.

When the Support is completely adjusted, except the tightening of the lacing in front, insert the pad (Fig. 2), as indicated in illustration, between the lowest rib and the hip-bone; then lace up and button the band fitted at the side, so as to cross the pad and keep it in position. In addition to this, the pad may be fastened with a safety-pin from the outside. Prices in all qualities up to 40 inches are:—

Qualities kept in Stock.	A Lavender Grey Jean.	E** White Horsehair (Summer).	J On Dr. Jaeger's System Natural Lasting.
	List Nos. 721	722	723

A Quality stocked in sizes 24 ins. to 36 ins.


** This quality is not kept in Stock, but will be delivered in six to seven days from receipt of order.

NOTE.—Instructions for measuring and extra charges for sizes over 40 inches are given on page 12.

INSTRUCTIONS FOR MEASUREMENT, ETC.

How to Measure:

The measurement should be taken over the skin, and, as indicated in the drawing, *round the line of greatest dilatation of abdomen (line 1-2), not round the waist*; this is quite sufficient.

**How to Order:**

If the measurement of the body consists of an uneven number of inches, deduct 3 inches; if of an even number, deduct 4 inches, to arrive at the proper size of Belt; because the ends of the Belt should not meet, but should allow sufficient room for expansion of the elastic and contraction of the body. Examples: If measurement be 33 ins., order a 30 ins. Belt. If measurement be 38 ins., order a 34 ins. Belt.

ALWAYS STATE WHEN ORDERING IF DEDUCTIONS HAVE BEEN MADE.

Stock Sizes:

All Belts catalogued are kept in sizes 28 ins. to 44 ins., unless otherwise stated.

Extra Charges:

All Belts described on pages 5 to 10, and all Belt Bandages and Body Belts made on Dr. Jaeger's System, when over 40 inches in size, are subject to an extra charge, according to the following scale:—

⁴²
1/-

⁴⁴
1/6

⁴⁶
2/-


⁴⁸
2/6

^{50 inches.}
3/-

Special Orders not returnable:

These, in addition to size of circumference, should carefully state depth in front as well as any abnormal requirements. From six to eight days are required for making Belts to special order, and an extra charge of 10 per cent. on the ordinary price is made. Belt Corsets to special order are also charged 10 per cent. extra when further measurements are stated beyond size of waist.

HOW THE BELT SUPPORT SHOULD BE WORN.



These Instructions are important, and should be carefully followed.

When the Belt is applied as shown in the illustration, the bands are equally drawn until the Belt fits well to the body. The portion of the band which is perforated and eyeleted is then fastened to the corresponding button or hook at the back. Care should be taken to draw the Belt as low down in front as possible, as this is where the support is most needed. When this is done, the elastic side-straps are fastened on the buttons at either side, so as to give support from below. Having completely adjusted the Belt, the main band is passed under the hook in front, and is buckled or otherwise fastened according to the style of Belt. **The lacing in front is in all cases tightened last, and should be loosened before putting on the Support.**


N.B.—The ends of the Belt at the back should not meet, but be about four inches apart, (for reasons see page 12). The Belt may be worn under or over the chemise, or over the drawers. In the two latter cases the chemise must sit loosely above the Belt, as otherwise the movements of the body may draw up the chemise and thus cause discomfort.

"DOMEN" ABDOMINAL BELT BANDAGES.

The following Bandages have been designed to supply the want of a well-constructed and comfortably sitting Belt for GENERAL use.

The selection of designs and the various styles of execution will suit every requirement and taste.


The Abdominal Bandages are somewhat lower in price than the Abdominal Belt Supports ; they are **not** intended for cases where **special complaints** have to be dealt with.



Design No. 21.

Grey Lasting, Elastic Bands and Watch-spring Steels as indicated, with Buttons and Eyelet Holes at back, to buckle in front.

List No. 724.



Design No. 26.

Grey Lasting, White Lining, Elastic Insertions at sides, Elastic Bands with Buttons and Eyelet Holes at back, to buckle in front, Watch-spring Steels, as in No. 21.

List No. 725.


Design No. 40B.

DR. JAEGER'S SYSTEM.

Natural colour, strong Stockinet (a very soft and pliable material, pervious to heat and perspiration), with Elastic Binding and Insertions over hips, Elastic Bands with Buttons and Eyelet Holes at back, to buckle in front, Watch-spring Steels as shown.

List No. 726.


For extra charges for sizes over 40 inches see page 3.



Design No. 45.

Cream colour Lasting, White Lining, Elastic Insertions at sides, to buckle at back, Watch-spring Steels as shown.


List No. 727.




Design No. 50A.

Cream Colour, strong Stockinet with Elastic Binding, Lasting Insertions at front and Back. Elastic Insertions at sides, and Elastic Bands with Buttons and Eyelet Holes at back, to buckle in front. Watch-spring Steels as shown.

List No. 728.





Design No. 61.

Woven Cotton Elastic, grey, with Elastic Binding, Elastic Bands with Buttons and Eyelet Holes at back, to buckle in front.

List No. 729.

Very strong, and well adapted for heavy exercise, riding, etc.


SPECIAL RIDING BELT FOR LADIES.

Design No. 65 $\frac{1}{2}$.

Made of Silk Elastic, shape as Design 61, to buckle in front.


This Belt is made to order only, List No. 730.

For Evening Wear.

Design No. 62 $\frac{1}{2}$.

Elastic Web with Silk (white, pink, cardinal, or blue), same as No. 61, but of shallow depth in front, for wear below corset, and to brace up abdomen.

List No. 731.



Design No. 55B.

White Horsehair Netting (a cool material, for wear in summer and in warm climates), trimmed in accordance with DR. JAEGER'S SYSTEM, Elastic Bands with Buttons and Eyelet Holes at back, to buckle in front.

List No. 732.

for all Belts up to 40 in. ; for extra charges see page 12.

THE "DOMEN" BELT-CORSET

is an actual combination of a Belt and a Corset, and affords the wearer all the advantages of wearing the two articles separately, without the inconvenience thereby entailed.

The Belt, as shown in the illustration, presents a firm front, but has elastic sides over the hips, which not only provide sufficient elasticity, but also ensure perfect adaptation to the figure.

The main bands of the Belt cross at the back in a similar manner to that described and illustrated on page 13, and being fastened by eyelet holes, are passed round the body, under the hooks in front, and are buckled to the required tightness.

Belt Corsets are not intended for special circumstances (such as are provided for on pages 6 to 10), but for general support, and for the improvement of the figure.

They are made in the following materials:—

List No.

733. A.—Grey Jean.

734. W.—White Coutile.

735. E.—White Horsehair.

736. D.—Black "Diamond Drill" Speciality.

737. J.—Nat. Color Pure Wool on Dr. Jaeger's System.

Special orders:—The measurements required for Belt Corsets made to special order, in addition to that described on page 12 (such Corsets are charged 10 per cent. extra, and cannot be returned to stock), are as follows:—

1. Waist Measurement.
2. Height under Arm from Waist to top of Corset.
3. Height at Back from Waist " "
4. Height in front from Waist " "
5. Chest Span.

Stock sizes are regulated by the waist measurement only.



"DOMEN" "STOOP-CURE."

For Illustration and Prices see next page.

THIS admirable appliance is a prompt and efficient corrective of tendency to stooping, and is so fitted to the figure beneath the dress as to be wholly unobservable.

It causes no discomfort to the wearer, who soon ceases to notice its presence, except when it gives a gentle reminder on any occasion of relaxing into the stooping habit.

Extraordinarily successful results have been achieved by the "DOMEN" "STOOP-CURE," and it possesses the immense advantage that it is readily and willingly worn, being neither unsightly nor uncomfortable, and, in fact, free from every objection whatsoever.

The "Stoop-Cure" is constructed of an anatomically shaped back-piece, containing watch-spring steels running parallel with the spine, and two springs extending from the centre line to the shoulders, with straps which draw the shoulders backwards, and thus increase the size of the thoracic cavity.

The "Stoop-Cure" adapts itself to every kind of weakness of the back, and to-day holds the field as the only perfect appliance to remedy round shoulders, and to maintain and preserve a correct and elegant carriage.

The "Stoop-Cure" is made in six sizes, each of which can be readily lengthened or widened to fit the wearer perfectly.

The only measurement required is that of the distance between the seams of the coat—or dress sleeves, across the shoulder blades; the mainspring should be appreciably smaller than such measurement.

"DOMEN" "SHOULDER-SHAPE."

THE "DOMEN" "SHOULDER-SHAPE" consists of two springs extending from the centre line of the back to the shoulders, with straps fastening under the arms and so preventing the shoulders from dropping forward. It has been merely designed to supply the want of a cheap appliance to keep the shoulders in the right position, but it is not intended for cases where stooping has already become a habit, or when support to the spine is necessary.

		SIZES	1.	2.	3.
Ordinary Quality	...	...	738	739	740
Extra	„ (Jaeger System)	...	741	742	743

Measurement same as for "Stoop-Cure."



PRICES.

Qual.	SIZES ...			0	1	2	3	4	5
L	Leather Back ...	...	List Nos.	744	745	746	747	748	749
S	White Horsehair do. ...	...	...	750	751	752	753	754	755
	On Dr. Jaeger's System	...	...	756	757	758	759	760	761

For Description see opposite page.


THE "DOMEN" MEN'S BELTS

have been designed to supply the want of a well-constructed and comfortably sitting belt for men's use. They not only give warmth to the abdominal regions, but also a support to the abdomen, which is much needed when riding, walking, or during any active exercise.

Design No. 403.

Natural colour, on Dr. Jaeger's system, with elastic insertions at sides, to buckle in front. List No. 738.

For extra charges for sizes over 40 inches see page 12.



Design No. 305.

Natural colour, on Dr. Jaeger's system, same shape as above, without elastic insertions at sides, to buckle in front. List No. 769.

For extra charges for sizes over 40 inches see page 12.

Design No. 32b.

Dr. Jaeger's System, Men's Girderbelt, natural colour pure Woollen Material, with Watch-spring Steels and Buckles to fasten at back. A most suitable light bandage for warmth and support. For extra charges for sizes over 40 inches see page 12. List No. 770.

Design 303.

A strong Felt Bandage, with Watch-spring Steels to buckle at back; made in grey or red colour. List No. 771.

Design 400.

A similar Belt, Pure Wool Felt, on Dr. Jaeger's System. List No. 772. For extra charge for sizes over 40 inches see page 12.

"DOMEN" CHOLERA BELTS.

For Protection against Chill and to maintain warmth over the Abdomen—not for Support.

N.B.—The following Belts are only made up to 40-in. in size.

No. 205.

Red Wool to fasten in front.

List No. 762.

(N.B.—The band on the right side passes through slot in the Belt at opposite end.)

No. 215.

Red Wool to buckle at back, with Japanned Steel Watch-spring Stays inserted at front, sides and back of the Bandage.

List No. 763.

No. 230.

Grey Flannel, with Buttons, and Eyelet Holes at back. The leather ends of the bands button to the side of the belt.

List No. 764.

"DOMEN" UMBILICAL BANDAGES.

FOR CHILDREN.

Measurement required:—Circumference of the Patient's body at the Navel.

In these bandages the elastic insertions, which easily become spoiled, are superseded by the use of japanned steel watch-springs, ensuring an agreeable and elastic adaptation of the bandage to the body of the wearer, preventing irritation and chafing even to the most sensitive skin, and, at the same time, making a displacement absolutely impossible. These bandages meet every requirement in the most simple and practical manner.

10-13.

14-16.

17-20.

List No.


...

...

765

766


767



I N cases of Varicose Veins the most effective remedy is a well-made elastic support, giving equal pressure to the limbs, and thus assisting the circulation of the blood. For this purpose the following anatomically-shaped Elastic Stockings, etc., have been designed. They are made of the best and most suitable web, can easily be washed in tepid water, and should be dried in the open air. They should not be worn during sleep, as muscular action must not be impeded.

HOW TO MEASURE.

The measurement should be taken over the skin, and in parts as indicated in drawing. As the limbs are liable to become swollen during the day it is best to measure in the morning. The exact measurement should be given, circumference and length, in inches; for stockings reaching over and below the knee, it is necessary to give the length above and below the knee separately, viz.:—For Stocking over knee give A-B-D-F.



STOCKING	...	...	A—D
Over Knee...	...	...	A—F
Over Thigh...	...	...	A—G
CALFPIECE	...	...	B—D
Over Knee...	...	...	B—F
Over Thigh...	...	...	B—G
THIGHPIECE	...	...	F—G
KNEECAP	...	...	D—F
Over Thigh...	...	...	D—G
ANKLET	...	...	A—B

Prices will be found in Separate List.

The ELASTIC STOCKINGS are made in the following qualities:—

- I.—Superior Fine Silk Elastic.
- II.—Medium Quality Silk Elastic.
- III.—Ordinary Quality Silk Elastic.
- IV.—Best Quality Cotton Elastic.

LADIES' SPECIALTY.

"DOMEN" OCCASIONAL BELT.

Reprinted from "Beauty's Queens."

"Yet another speciality is the 'OCCASIONAL BELT,' which we cannot too highly estimate for ladies' occasional wear, possessing, as it does, so many advantages over the many primitive and inefficient appliances in use. It is so constructed as to entirely remove all possible source of discomfort and risks of disease and chill, its completeness being of special advantage for travelling. Without describing its component parts, we would merely advocate its adoption, and emphasize the fact that many of the evils from which many delicate women suffer may be avoided by its use."

List No. 778-787.



ASK FOR SPECIAL LEAFLETS.





"DOMEN"
REG^d